



Operating Light ZOL-A30

Operating Light ZOL-A30

Operating light ZOL-A30 features a high quality reflector allowing the bulb to focus concentrate the light to its specific operational area. Designed to use by medical personnel for individual medical examination and small procedures. Failure of the primary light will be automatically replaced by standby light.

Features:

High quality reflector

Automatic switch over between the primary and standby lights

Brightness focus (light diameter density) can be customized

Applications:

Used in hospitals to conduct modern surgical operations as ideal light source.

Specifications :

Illumination Intensity	50,000 lux ~ 100,000 lux
Lamp Holder Diameter	500 mm
Light Field	100 mm ~ 280 mm
Temperature Increased	? 12 °C
Depth of Illumination	? 700 mm
Bulb Voltage	24 V (AC)
Mains Voltage	220 V ± 20 % <hr/> 50 Hz
Temperature Increased of Doctor's head	? 2 °C
Power Input	220 VA



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA
Email: info@zimed.com | Website: Zimed | Phone: +1 302 786 1942