



**Ventilator machine ZVM-A40**

## Ventilator machine ZVM-A40

Ventilator machine ZVM-A40 is a multipurpose ventilator used for emergency clinical resuscitation and respiratory therapy for adult and neonate patients. There is backup power supply facility during power failure and Alarm system give alert sign in case of any failure in the functions.

### Features:

Multipurpose : adult and neonate

Ventilation mode-1) Adult: IPPV(VC)/(PC), SIMV, SIMV+PSV, CPAP, PSV 2) Neonate: IPPV (PC), IMV, CPAP continues flow, limited pressure, Time cycle 3) Emergency: IPPA (PC) constant pressure

Inspiratory flow rate: 1) Adult: 10 to 120 L/min (inspiratory phase) 2) Infant: 4 to 24 L/min (continuous)

Respiratory rate, minute volume, tidal volume are visible in the ventilator monitor

Maximum, mean and plateau pressures are observed in the airway pressure monitor

Presence of alarm system

Backup power supply

### Applications:

Used in emergency care units and also in intensive care units of hospitals.

## Specifications :

<b>Respiratory rate</b>	IPPV: 5 to 60 bpm, SIMV: 1 ? 12 bpm
<b>Inspiratory time</b>	0.2 ? 3 sec
<b>Minute volume</b>	0-40 L/min
<b>Pressure limit</b>	1 to 9 kPa
<b>Adjustable oxygen concentration</b>	21-100 %



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA  
Email: [info@zimed.com](mailto:info@zimed.com) | Website: [Zimed](http://Zimed.com) | Phone: +1 302 786 1942