



Electrosurgical unit ZEU-A12

Overview

Electrosurgical unit ZEU-A12 is used for bipolar surgeries. It offers energy output up to 150 W. Therapeutic methods includes bipolar coagulation with lower thermal damage risks. It offers separate control panel for a range of operations.

Features:

- Energy output up to 150 W
- Digital display for visualization of working parameters
- Push button for quick change of parameters
- Error alarm for failure detection system
- Easy to operate

Applications:

This is used during general surgery. Also used during neurosurgery.

Specifications :

Output Power	80 W
Power Rating	300 VA \pm 10 %
Operating frequency(monopolar)	330 kHz/430 kHz
Operating frequency(bipolar)	1000 kHz
Bipolar coagulation	0 to 80 W (300 ohm)
Biopolar	0 to 80 W (100 ohm)



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA
Email: info@zimed.com | Website: www.zimed.com | Phone: +1 302 786 1942