



Manual Orthopedic Traction Bed ZOB-A101

Overview

Manual orthopedic traction bed ZOB-A101 is designed with a strong, fireproof ABS panel with 360° rotation. Solid bed base made of 5.0 mm stamped rotating plate and retractable hand crank are also equipped in the bed. It consists of TPR wire with an anti-winding hard shell, an L-shaped folded guard rail, and a bended tube for the safety of patients. The bed is equipped with an IV pole and drainage hooks as well.

Features:

Safe, simple, reliable, and easy-to-use

IV pole and drainage hooks

Solid bed base

ABS strong fire-proof panel

Head and foot board

5” hard-shell castors with break

Traction frame

Applications:

Orthopedic beds are used for patients in the fields of hospitals, medical institutes, etc.

Specifications :

Trendelenburg/ Reverse	0 - 12°
Knee – rest adjustment	0 - 35° (±10°)
Back – rest adjustment	0 - 75° (±10°)
Mattress platform	1925 × 830 mm
Bed height adjustment	460 – 710 mm
Castor	4 pcs
IV pole prevision	4 pcs
Drainage Hook	2 pcs
Traction Frame	1 pc
Manual crank	4 pcs
Dimensions (L × W)	2180 × 960 mm
Safe Working Load	250 Kg



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA
Email: info@zimed.com | Website: www.zimed.com | Phone: +1 302 786 1942