



## **Electric Orthopedic Traction Bed ZOB-A201**

## **Electric Orthopedic Traction Bed ZOB-A201**

Electric orthopedic traction bed ZOB-A201 consists of a stamped-form platform with no welding scar. Environmental plastics can also be sterilized at 100°C and have a tensile strength of 30 MPa. It features anti-bacterial steel with a paint thickness of 0.12 mm and a brightness of 60°. Double-sided castors with a central locking system, a hand controller, and a linear motor are incorporated in the system. The bed is equipped with ABS side rails with extra height and a safe lock.

### **Features:**

It is simple, safe, reliable, easy-to-use  
Hand controller, angle indicator  
ABS side rails and bed ends with extra height and safe lock  
Double side castors with central locking system and pedal  
IV pole provision, drainage hook  
Trendelenburg, height adjustment  
Back-rest and knee-rest adjustment too  
Can bear up to 250 kg

### **Applications:**

Orthopedic beds are used for patients in the fields of hospitals and medical institutes.

## Specifications :

<b>Back-rest adjustment</b>	0 - 75° ( $\pm$ 10°)
<b>Knee-rest adjustment</b>	0 - 35° ( $\pm$ 10°)
<b>Height adjustment of use base</b>	435 – 725 mm
<b>Trendelenburg</b>	0 - 12°
<b>Reverse Trendelenburg</b>	0 - 12°
<b>Double Side Castors</b>	4 pcs (central locking system)
<b>Mattress Platform</b>	1935 × 900 mm
<b>IV pole prevision</b>	4 pcs
<b>Drainage Hook</b>	2 pcs
<b>Dimensions (L × W)</b>	2250 × 1060 mm
<b>Working Load</b>	250 Kg



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA  
Email: [info@zimed.com](mailto:info@zimed.com) | Website: [www.zimed.com](http://www.zimed.com) | Phone: +1 302 786 1942