



## **Infant Incubator ZII-A13**

## **Infant Incubator ZII-A13**

Infant Incubator ZII-A13 is employed with the latest technology to provide a secure environment in which infants can develop as their organs with observation and care. Uses temperature range from 20° to 30° with humidity range from 30% to 75%. Its transparent cover allows clear observation of the infant without sacrificing the controlled environment. Port holes are sealed with rubber gloves allowing the nurse to carry a baby without contamination.

### **Features:**

- Uses two Control modes: Air mode and Skin mode
- Equipped with a large display and audible & visual alarm
- Provided with four tubing grommets for the safety
- Can warm up the cabinet within 45 minutes
- Includes four poles, protective guard rail, and storage cabinet
- Allows keypad lock, and self-test function
- Integrated with an alarm that includes an alarm silence function to prevent protection against power, fan, sensor, and system failure
- Bears comfortable mattress and hood with dual slope and 6 hand-ports

### **Applications:**

Infant incubators are widely used for maintaining the desired temperature for premature babies in the neonatal intensive care unit (NICU), and pediatric department.

## Specifications :

<b>Working Temperature</b>	20? to 30?
<b>Working Humidity</b>	30% to 75%
<b>Air Velocity</b>	0.3 m/s
<b>Air Temperature Control Range</b>	20.0? to 37.0?
<b>Skin Temperature Control Range</b>	20.0? to 37.0?
<b>Noise Level</b>	? 55 dB
<b>Temperature Fluctuations</b>	? 0.5?
<b>Operation Port</b>	6
<b>Front Access Panel</b>	826 × 310 mm
<b>Mattress Tray Size</b>	630 × 360 mm
<b>High Air Temperature Alarm</b>	+3.0?
<b>Low Air Temperature Alarm</b>	-3.0?
<b>High Skin Temperature Alarm</b>	+1.0?
<b>Low Skin Temperature Alarm</b>	-1.0?
<b>Air Temperature Alarm</b>	?38.0?; ?40.0?
<b>Skin Temperature Alarm</b>	?38.0?; ?39.0?
<b>Dimensions (L×W×H)</b>	116 ×56 ×137 cm
<b>Weight</b>	87 kgs



<span>Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA</span>  
Email: [info@zimed.com](mailto:info@zimed.com) | Website: [www.zimed.com](http://www.zimed.com) | Phone: +1 302 786 1942