



Ophthalmoscope ZOS-A10

Overview

Ophthalmoscope ZOS-A10 features 19 adjustable corrective lenses, spanning a comprehensive range from -20D to +20D. It is equipped with a deluxe LCD screen handle that allows users to adjust the brightness of the light source. The built-in current IC circuit ensures a stable light effect, providing a warm white illumination. Our Ophthalmoscope's extensive range of lenses ensures precise and customizable vision correction for a variety of diagnostic needs.

Features:

- Optical glass lens
- 5 different apertures
- Optical SMT LED technology
- Low power consumption

Applications:

Our Ophthalmoscope is used in clinical settings, including ophthalmology and optometry, for examining the retina and diagnosing a range of eye conditions.

Specifications :

Illumination	Direct illumination
Magnification	15x
Corrective Lens	19, adjustable
Range of Lenses	-20D to +20D
Bulb	LED 3.0Vdc 80mA
Bulb Model	206L
Luminance	14 lumens
Color Temperature	3200K \pm 5% for warm white light
Operation Temperature	10? to 35?
Operation Humidity	30% to 75%
Operation Air Pressure	700hPa to 1060hPa
Storage and Transport Temperature	-20? to 70?
Storage and Transport Humidity	10% to 90%
Storage and Transport Air Pressure	500hPa to 1060hPa
Battery	2 \times AA alkaline batteries
Dimension	190 \times 110 \times 40 mm
Packaging Dimension	400 \times 240 \times 230 mm
Net Weight	6 kg
Gross Weight	6.5 kg



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA
Email: info@zimed.com | Website: www.zimed.com | Phone: +1 302 786 1942