



Patient Monitor ZPM-A21

Patient Monitor ZPM-A21

Patient monitor ZPM-A21 offer reliable monitoring technologies with accurate reading and optimized performance. Easy monitoring mode for accurate prediction of patient's condition. Auto save mode for temperature, ECG, respiratory and pulse rate.

Features:

Touch screen display

Display's temperature, ECG, respiration and pulse rate, SpO₂, NIBP

Compact and portable

Retains heart rate, body temperature, SpO₂, respiration rate status for 72 hrs.

Audible/ visual alarm system

Suitable for adult, pediatric and neonate patients

Real time battery usage

Optional built-in Printer, ETCO₂

Applications:

Used to monitor the parameters of ICU patient who need daily care.

Specifications :

Display type	LCD (TFT)
Display size	12.1 inch
Standard parameters	Temp, ECG, pulse, respiratory rate, SPO2, NIBP
Dimensions	430 × 300 × 400 mm
Weight	7 kg



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA
Email: info@zimed.com | Website: www.zimed.com | Phone: +1 302 786 1942