



Surgical Suction Pump ZSP-A20

Surgical Suction Pump ZSP-A20

Surgical suction pump ZSP-A20 features oil free pump, operate with flow rate more than 80 L/min helps to achieve increased performance and faster suction. Four containers are incorporated to collect fluid, air filter is added to prevent damage on pump. High capacity, polycarbonate jar is used. Sensitive operation with less noise.

Features:

- Four containers incorporated
- Flow rate can be above 80 L/min
- Overflow protection prevents fluid entering in to the inner pipe
- High intensity polycarbonate plastic jar
- Oil free piston pump

Applications:

Used in clinics, hospital where general surgeries conducted, Liposuction, endoscopy, wound drainage

Specifications :

| | |
|-------------------------------|-----------------------------------------|
| Max Vacuum Pressure | 680 mmHg |
| Vacuum Range | 0.02 MPa ~ 0.09MPa (680 mmHg) |
| Flow Rate | ? 80 L/min |
| No. of Container (Jar) | 4 |
| Noise | ? 60 dB |
| Input Power | 280 VA |
| Container Capacity | (4000 ml × 2) ----- (2000 ml × 2) |
| Power Supply | AC 220 V ----- 50 Hz |



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA
Email: info@zimed.com | Website: www.zimed.com | Phone: +1 302 786 1942