



**Ventilator machine ZVM-A50**

## **Ventilator machine ZVM-A50**

Ventilator machine ZVM-A50 is pneumatically driven and electrically controlled device. Tidal volume delivery ranges from 10 to 300 ml. It's newly designed to meet the special need of clinical respiratory management and respiratory therapy for infant and children with respiratory disorders.

### **Features:**

- Compact and highly reliable
- Innovative functions, convenient operation with low tidal volume
- Essential instrument in emergency and ICU rooms
- Audible/ visual alarm systems
- Backup power supply

### **Applications:**

Specially designed for new born babies (infant), used in clinical and intensive care units.

## Specifications :

<b>Tidal volume</b>	10 - 300 ml
<b>Ventilation mode</b>	CMV, IMV, CPAP, Manual
<b>Respiratory rate</b>	1 to 120 bpm
<b>I:E ratio</b>	1:0.1 - 1:9.9
<b>Minute volume</b>	? 5 L/min
<b>Adjustable oxygen concentration</b>	21 - 100 %



Zimed Healthcare Inc. 1400 Broadfield Blvd, Houston, Texas 77084 USA  
Email: [info@zimed.com](mailto:info@zimed.com) | Website: [www.zimed.com](http://www.zimed.com) | Phone: +1 302 786 1942